



A. Model A9024



B. Model A9087

## Common Features:

### TREATMENT TABLES –

MODELS A9002, A9024, A9041, A9087;

- 30" x 78" upholstered top with 2" High-Density foam.
- High-pressure laminate legs, apron and H-Brace.
- New enhanced leg hardware system.
- Legs have protective black vinyl scuff caps.

**COLORS:** (FOLD OUT FRONT COVER)

**Laminate** – Natural Oak or Wild Cherry

**Upholstery** – 9 Pro-Form vinyl colors

## ENHANCED LEG HARDWARE SYSTEM

Attach or detach legs quickly and easily with heavy-duty  $\frac{5}{16}$ " bolts into embedded metal threaded inserts for durable connection. Easily tighten and re-tighten as needed.



## A. Models A9002, A9024

### Treatment Table

#### SPECIAL FEATURES:

- A9024 has full width laminate shelf.

MODEL	LENGTH	WIDTH	HEIGHT
A9024	78"	30"	31"
A9002 (not shown)	same as A9024 w/o shelf		

## OPTIONS

MODELS A9002, A9024, A9041, A9073

\*08 – Drawer (for A9024, A9041 only)

\*69 – Gas Spring Backrest 0° to 75°

A946 – (1) Lock/Bolt for Double Door

\*FACTORY INSTALLED

SEE PAGE 4

## C. Model A9041

### Treatment Table

#### SPECIAL FEATURES:

- Laminate shelf and open storage cabinet.

MODEL	LENGTH	WIDTH	HEIGHT
A9041	78"	30"	31"

## B. Model A9087

### 600 lb. Trainers Table

#### SPECIAL FEATURES:

- (6) Legs and plywood top for added strength.
- Laminate shelf with (2) long stretchers.

MODEL	LENGTH	WIDTH	HEIGHT
A9087	78"	30"	31"

## D. Model A9073

### Cabinet Table

#### FEATURES:

- (3) Two-door cabinets with laminate interiors and adjustable shelf.
- Optional A946 Lock/Bolt for double door (specify).
- Natural Oak or Wild Cherry laminate.
- 9 ProForm vinyl colors (specify).

MODEL	LENGTH	WIDTH	HEIGHT
A9073	78"	30"	31"



C. Model A9041  
Shown with option:  
08 – Drawer



D. Model A9073