


A. MODEL A9068

Shown with options:
 69 - Gas Spring Backrest
 (3) A946 Lock/Bolts
 (1) A905 Embossed Logo



COMMON FEATURES

SPLIT LEG TABLES - MODELS A9064, A9065, A9068:

- Extremely versatile and functional split leg design.   
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.
- Comfortable 2" thick high-density foam upholstered tops.


COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak, Wild Cherry or Folkstone Gray

Upholstery - 10 Pro-Form™ vinyl colors

OPTIONS

MODELS A9064, A9065, A9068

*69 - Gas Spring Backrest 0° to 75° 

*99-L - Left hand table (A9068 only)

*A946 - (1) Lock/Bolt for Double Doors

*FACTORY INSTALLED

B. MODEL A9064



**A. MODEL A9068
SPLIT LEG TABLE**

SPECIAL FEATURES:

- Six hinged doors with interior shelf offer concealed supply storage.
- Solid finished back.

MODEL	LENGTH	WIDTH	HEIGHT
A9068	78"	30"	31"

**B. MODEL A9064
SPLIT LEG TABLE**

SPECIAL FEATURES:

- Six hinged doors with interior shelf on both sides of table.

MODEL	LENGTH	WIDTH	HEIGHT
A9064	78"	30"	31"

C. MODEL A9065



**C. MODEL A9065
SPLIT LEG TABLE**

SPECIAL FEATURES:

- Open shelf storage, accessible from both sides of the table.

MODEL	LENGTH	WIDTH	HEIGHT
A9065	78"	30"	31"