



A. MODEL A9069
Shown with option:
69 - Gas Spring Backrest

A. MODEL A9069
KNEE AND HIP FLEXION SPLIT LEG SIX LEG TABLE



- High-pressure laminate legs, apron, and H-Brace.
- Laminate shelf with (2) long stretchers.
- New enhanced hardware system
- Legs have protective black vinyl scuff caps.

MODEL	LENGTH	WIDTH	HEIGHT
A9069	78"	30"	32"



Straight leg lift

COMMON FEATURES

KNEE AND HIP FLEXION SPLIT LEG TABLES - MODELS A9069, A9070, A9071:

- Extremely versatile and functional split leg design.
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.
- Each padded leg rest can be adjusted with hip and knee joints locked into 45° angle with easy to operate hinge mechanism.
- Full plywood subtop for added strength.
- 30" x 78" upholstered top with comfortable 2" thick high-density foam.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak, Wild Cherry or Folkstone Gray (Folkstone - A9070, A9071 only)

Upholstery - 10 Pro-Form™ vinyl colors

B. MODEL A9070

Shown with option:
69 - Gas Spring Backrest



B. MODEL A9070
KNEE AND HIP FLEXION SPLIT LEG OPEN STORAGE TABLE



- Open shelf storage, accessible from both sides of the table.
- Full plywood subtop for added strength.

MODEL	LENGTH	WIDTH	HEIGHT
A9070	78"	30"	32"

C. MODEL A9071

Shown with option:
69 - Gas Spring Backrest



C. MODEL A9071
KNEE AND HIP FLEXION SPLIT LEG CABINET TABLE



- Six hinged doors with interior shelf offer concealed supply storage.
- Full subtop for increased strength.

MODEL	LENGTH	WIDTH	HEIGHT
A9071	78"	30"	32"

OPTIONS

MODELS A9069, A9070, A9071

*69 - Gas Spring Backrest 0° to 75°

*99-L - Left hand table (A9071 only)

*A946 - (1) Lock/Bolt for Double Doors

*FACTORY INSTALLED