



COMMON FEATURES

S&W TREATMENT TABLES:

- Natural wood finish is standard on all S&W Tables.
- H-Brace Tables have triple-bolted, solid hardwood legs.
- Cabinet tables are quality hardwood plywood and veneer construction.
- Extra-firm 2" high-density foam tops with rounded corners for comfort.









COLORS: (SEE INSIDE FRONT COVER)

Upholstery – 10 Pro-Form™ vinyl colors



MODEL ATT2472


MODEL	LENGTH	WIDTH	HEIGHT
ATT2472 	72"	24"	31"
ATT2478 	78"	24"	31"
ATT2772 	72"	27"	31"
ATT2778 	78"	27"	31"
ATT3072 	72"	30"	31"
ATT3078 	78"	30"	31"

OPTIONS

AFW – Finished Wood Shelf

AU – Upholstered Shelf

27 – Paper Dispenser/Cutter Combo

***69** – Gas Spring Backrest 0° to 75°
(N/A on 24"W Tables) 

***A946** – (1) Lock/Bolt for Double Doors

***FACTORY INSTALLED**



MODEL ATT2772

Shown with option:
AU – Upholstered Shelf

A. MODEL ATT314

SPECIAL FEATURES:

- Four door cabinet with center divider and one adjustable shelf on each side.
- Leg levelers for uneven floors.
- Satin finish wire pulls.



MODEL	LENGTH	WIDTH	HEIGHT
ATT314	72"	30"	31"

B. MODELS ATT385, ATT386

SPECIAL FEATURES:

- Includes 30" long Gas Spring Backrest, 0° to 75°, as standard feature.
- Includes full-width AU upholstered shelf.



MODEL	LENGTH	WIDTH	HEIGHT
ATT385	72"	30"	31"
ATT386 (not shown)	78"	30"	31"



AFW Wood Shelf

AU Upholstered Shelf



A. MODEL ATT314

Shown with option:
69 – Gas Spring Backrest



B. MODEL ATT385

Shown with standard
69 – Gas Spring Backrest