

# INVACARE® SHOWER CHAIRS



---

**Invacare® I-Fit™** and **Invacare® CareGuard™ Shower Chairs** offer stability, comfort and tool-less assembly in seconds. Invacare® shower chairs are designed for individuals who require seating support when bathing or showering and provide comfort and safety in the bathtub. Available with a composite or vinyl padded seat, shower chairs fit in most standard tubs and showers.

---

---

### Invacare® I-Fit™ Shower Chair Features

Model # 9781, 9780

Tool-less assembly

Deep and wide seating surface


Contoured seat with drain holes and ergonomic backrest

Plastic legs with twist height adjustments

Numbering system built into legs to help keep chair level

---

#### HCPCS Code: E0240



1. Tool-less assembly

2. Retail packaging



9781: With Back



9780: Without Back

---

### Invacare® Vinyl Padded Shower Chair Features

Model # 9872


Vinyl seat and back are padded with 1.25" open-cell foam for comfort

Nine height adjustments in half-inch increments

Non-marring rubber tips

---

#### HCPCS Code: E0240



1. Non-marring rubber tips



9872

---

### Invacare® CareGuard™ Tool-Less Shower Chair Features

Model # 96-2, 91-2

Tool-less assembly

Deep and wide seating surface

Contoured seat with ergonomic backrest

Textured seat with drain holes



Hand hole built into seat back for easy transport

Aluminum legs with half-inch snap button adjustments to help with stability

Wide-based crutch tips to help increase stability

---

#### HCPCS Code: E0240



- 1. Tool-less assembly
- 2. Retail packaging



### Invacare® Bariatric Shower Chair Features

Model # 9785-2

700 lb. weight capacity

Wide, deep seating surface

Contoured seat and backrest

Adjustable height legs

Durable aluminum frame with dual cross brace for added stability

Backrest easily snaps into assembled seat frame with no tools required

---

#### HCPCS Code: E0240



- 1. Bariatric size



**Invacare Corporation**

www.invacare.com

**USA**

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

**Canada**

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
®. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Specifications are subject  
to change without  
notification.  
Form No. 11-237  
111267

**Yes, you can.®**

Model #	Overall Dimensions	Overall Base Dimensions	Seat to Back Height	Product Weight	Weight Capacity	Units Per Carton	Limited Warranty
9780	18" D 20" W 15"-21" H	17.25" D 16.25" W	N/A	8 lb. 4 oz.	400 lb.	2	3 Years
9781	18" D 20" W 19.25"-35.25" H	17.25" D 16.25" W	14.25"	10 lb.	400 lb.	2	3 Years
9872	16" D 15.5" W 17.5"-22" H	16.25"-17" D 19.25"-19.75" W	14.75"	9 lb.	250 lb.	2	3 Years
91-2	12.75" D 21" W 16.125"-22" H	9.5"-18.875" D 17.25" W	N/A	5 lb. 14 oz.	315 lb.	2	3 Years
96-2	12.75" D 21" W 16.125"-22" H	13.75"-18.875" D 17.25" W	17.25"	6 lb. 6 oz.	315 lb.	2	3 Years
9785-2	16.25" D 17" W 16.25"-20" H	21"-22" D 18.25"-18.75" W	13.5"	11 lb.	700 lb.	2	3 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.