

INVACARE®  
**HOMEFILL®**  
CYLINDERS



Invacare offers a wide variety of **Invacare® HomeFill® Cylinders** to fit your client's lifestyle and daily oxygen needs. These cylinders range in size from as small as a water bottle to large enough to last almost 16 hours\*.

Now HomeFill cylinders are also available with CGA870 compatible post valves and integrated conservers. With a full line of cylinders, the Invacare HomeFill Oxygen System meets the needs of most oxygen patients.

\*Duration Varies Depending on the User's Oxygen Needs.

## Integrated Conservor Cylinders

HCPCS Code: K0738



HF2PCM2KIT, HF2PCL4, HF2PCL6, HF2PC9



HF2ICL4, HF2ICL6, HF2IC9,

## Post Valve Cylinders

HCPCS Code: K0738



HF2POST6



HF2POST9



HF2POSTD

## Continuous Flow Cylinders

HCPCS Code: K0738



IOH2MD

### Conserved Duration (hrs. @ 20 BPM)

Setting	1	2	3	4	5	6
HF2PCM2KIT	3.4	2.2	1.6	1.3	1.1	N/A
HF2PCL4	6.7	4.3	3.2	2.6	2.1	N/A
HF2PCL6	9.8	6.3	4.7	3.8	3.1	N/A
HF2PC9	15.2	10.5	7.3	5.9	4.8	N/A
HF2ICL4	5.9	2.9	2.0	1.5	1.2	1.0
HF2ICL6	8.9	4.5	3.0	2.2	1.8	1.5
HF2IC9	12.8	6.4	4.3	3.2	2.6	2.1

## Invacare Corporation

www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2013

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
©. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Specifications are subject  
to change without  
notification.  
Form No. 07-066  
Rev. 06/13  
130343



Yes, you can.®



Model #	Height	Diameter	Weight	Continuous Flow Setting	Capacity	Duration @ 2 LPM Flow	Fill Time
HF2PCM2KIT	8.25"	3"	2 lb.	2 LPM	57 L	24 min.	24 min.
HF2PCL4	9"	4.38"	3.1 lb.	2 LPM	113 L	55 min.	55 min.
HF2PCL6	11"	4.38"	3.5 lb.	2 LPM	170 L	75 min.	75 min.
HF2PC9	14"	4.38"	4.8 lb.	2 LPM	255 L	125 min.	125 min.
HF2ICL4	9"	4.38"	3.2 lb.	2 LPM	113 L	55 min.	55 min.
HF2ICL6	11"	4.38"	3.6 lb.	2 LPM	170 L	75 min.	75 min.
HF2IC9	14"	4.38"	4.9 lb.	2 LPM	255 L	125 min.	125 min.
HF2POST6	12"	4.38"	4.2 lb.	N/A	170 L	75 min.	75 min.
HF2POST9	15"	4.38"	5.4 lb.	N/A	255 L	125 min.	125 min.
HF2POSTD	21"	4.38"	6.5 lb.	N/A	422 L	3.5 hrs.	3.5 hrs.
IOH2MD	19.75"	4.38"	7 lb.	0-6 LPM	422 L	3.5 hrs.	3.5 hrs.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.