

INVACARE®  
**JASMINE™**  
PREMIER SERIES FULL BODY LIFT



The **Invacare® Jasmine™ Full Body Lift** has optimum functionality and innovative features that provide safer patient handling within a 500 lb. weight capacity. Caregivers can rely on safer operation and patients on maximum comfort during everyday transfers. Beyond that, the Jasmine lift expertly functions in more challenging situations with a rotating 360° 6-point hanger bar and 16.4" minimum height range, essential for assisting patients who may have fallen to the floor.

The stable spring-loaded hanger bar produces an audible click when the sling is engaged. This assures the caregiver that the sling is properly attached. That means less stress for the caregiver, a stable transfer for the patient and less likelihood of injury for both. Optimized maneuverability is key for a comfortable transfer. Jasmine lifts feature a battery-powered boom and wide-opening base, so caregivers can easily walk behind them and wheelchairs can be accommodated up front. The push handle accommodates any caregiver, and with a lighter-weight aluminum and steel construction, control in movement is straightforward and trouble-free.

## Features

**Lifting Range:** This lift expertly functions in more challenging situations with a rotating 360° 6-point hanger bar and 16.4" minimum height range, essential for assisting patients who may have fallen to the floor

**Maneuverability:** The push handle accommodates any caregiver. And with a lighter-weight aluminum and steel construction, control in movement is straightforward and trouble-free

**Ease For Caregivers:** Jasmine lifts feature a battery-powered boom and wide-opening base, both controlled through the hand pendant

**Security:** The stable spring-loaded hanger bar produces an audible tone when the sling is engaged. This assures the caregiver that the sling is properly attached

Cradle

XL Cradle

Scale : JLS5

Anti-Entrapment motor safety device



**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
©. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Specifications are subject  
to change without  
notification.  
Form No. 11-243  
111279



1. 360° 6 Point Hanger Bar
2. Hand Pendant
3. Push Handles
4. Cradle
5. Scale



Model #	Height At Sling Hook Up	Base Dimensions	Casters	Weight	Weight Capacity	Low Battery Alarm	Emergency Controls
JASMINE	16.4"-74.6" H	4.8"H Clearance 46.7"W Open 25.8"W Closed 49"L	4"	98 lb.	500 lb.	Audio	Stop Up and Down
	Battery	Charge Time	Charger Input	Charger Output	Limited Warranty		
	24V DC Rechargeable 2 w/Unit	6 hr.	100-240V AC	29.5 DC	3 Years Lift 2 Years Electronics		