



Hydraulic Patient Lifts

LIFT, PATIENT, HYDRAULIC, 6 POINT CRADLE

- Elevating hydraulic lift.
- Base opens with an easy hand lever.
- Lift comes with a full-body sling packaged inside box.
- 6-point cradle type.

Packaging	1 Each / Each
Manuf / Supplier	Medline
Manuf / Supplier #	MDS88200D

SKU MDS88200D

Specifications

Specifications	Values
Latex Free	Yes
Product Type	Hydraulic Lift
Weight Capacity	400 lb

Packaging Details

Unit of Measure	Conversion	Net / Gross Weight (lbs)	Volume (cubic ft)	Dimensions (inch) L x W x H	GTIN
Each (EA)	1.0 Each	0.0 / 93.156	8.865	46.0 x 12.0 x 27.75	10080196764259

HCPCS Information

HCPCS Code: E0630

HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.

Generated 08/31/2012