

## Truly interoperable connectivity

Mindray's flexible connectivity options enable you to future-proof your hospitals systems, for scalable growth and future innovation. Our vendor neutral and comprehensive connectivity solutions make digital patient data integration simple and secure.

- Seamlessly connects to your EMR systems
- Enabled caregiver login and patient identification
- Hybrid BeneVision central monitoring system to maximize your investment
- M-IoT device manager provided to screen the clinical efficiency of the equipment



## VS 9

Vital Signs Monitor

## True to Life

## True to form quality and reliability

Created with longevity in mind, the VS 9 is supported by innovative maintenance software and a robust design.

### Trouble-free Cleaning

Validated with up to 51 leading cleaning agents, to aid faster wipe-downs and less degradation over time.

### Easier installation

Step-by-step configured for your own clinical requirements with a smartphone-style Startup Wizard.



Useful workboard with optional extended battery for longer runtime

Tool-free design for faster installation under one minute

Dedicated accessory management as you prefer

Quiet and smooth casters for ease of use





## See the true picture of patient deterioration

To help ensure a timely reaction to patient deterioration and greater compliance with hospital guidelines, the VS 9 provides a comprehensive range of vital signs parameters with a variety of tools to support clinical decisions.

**Automatic EWS**

- Fully compliant with all EWS standards
- Provide customizable escalation messages

**Orthostatic BP**

- Innovatively help to assess and avoid fall risk
- Customizable due to your own clinical requirements

**BP Average**

- Practice unattended AOBP by your own workflow
- Provide accurate values for real patient status

From inpatient admissions to spot checks, the VS 9 Vital Signs Monitor unlocks the full potential of routine observations with comprehensive and flexible workflow options empowering clinicians to deliver high quality bedside care to each and every patient.


## True-to-life vital signs

The VS 9 delivers precision vital signs measurements at rapid speeds, maximizing patient comfort and clinical efficiency.

Mindray TrueBP™ inflation NIBP measurement in 15 seconds

SpO<sub>2</sub> measurement options of Mindray, Masimo SET and Nellcor

ARTEMA Sidestream CO<sub>2</sub> with 50ml/min sampling rate



## Truly flexible workflows

Providing comprehensive documentation at the point of care, the VS 9's workflows unlock the full potential of bedside observations and help consolidate and speed up bedside assessments.

### Vivid user interface

- 10.1" high-resolution screen to view more information
- Intuitive gesture operation with capacitive touchscreen

### Capture more at the bedside

- Up to 30 customizable manual inputs
- Infographic pain charting tool

### Create the right workflow

- Customizable screen layout and parameter settings
- Configurable workflows by department profiles
- Improve caregiver compliance to hospital guidelines

