

3M™ Red Dot™ Pre-wired Monitoring Electrodes



# Wired for Reliability

## An option for reducing HAIs

Our 3M™ Red Dot™ Pre-wired Monitoring Electrodes, available in adult, neonatal and infant sizes, are designed with both clinician convenience and patient comfort in mind. As a single-use disposable product, they can potentially reduce the risk of Healthcare Associated Infections (HAIs) associated with cross-contamination caused by dirty or insufficiently cleaned reusable leadwires.

## Traveling between departments is no longer a disconnect

Used in conjunction with the YokeMate LWS® cable and adapter system, one set of Red Dot™ Pre-wired Monitoring Electrodes can travel with your patients between departments throughout their hospital stay. No matter what ECG monitors are being used, the YokeMate LWS® system requires adapters only 10% of the time – much less than other systems.



Contact 3M for help in choosing appropriate cables and adapters to complement the monitors used in your facility.

3M Health Care Helpline 1-800-228-3957  
[www.3M.com/reddotelectrodes](http://www.3M.com/reddotelectrodes)



## Ordering Information

### Adult Pre-wired Monitoring Electrodes

Product No.	Configuration	Leadwire Length	Gel Type	Packaging
2289PA	5-leadwire	28" leadwire	Conductive Adhesive	500 electrodes/case, 5/bag, 100 bags/case
2289PAL	5-leadwire	48" leadwire	Conductive Adhesive	500 electrodes/case, 5/bag, 100 bags/case
2289T	3-leadwire	28" leadwire	Conductive Adhesive	300 electrodes/case, 3/bag, 100 bags/case

### Neonatal and Infant Pre-wired Monitoring Electrodes

Product No.	Configuration	Size	Leadwire Length	Backing Type	Gel Type	Radiolucent	Packaging
2269T	3-leadwire	1.25" x .75"	22" leadwire	Clear Tape	Conductive Adhesive	Yes	300 electrodes/case, 3/strip, 10 bag/box, 10 box/case
2282	3-leadwire	.9" x .9"	18" leadwire	Soft Cloth	Conductive Adhesive	Yes	300 electrodes/case, 3/bag, 100 bag/case
2282L	3-leadwire	.9" x .9"	26" leadwire	Soft Cloth	Conductive Adhesive	Yes	300 electrodes/case, 3/bag, 100 bag/case
2283	3-leadwire	1.3" x .9"	18" leadwire	Soft Cloth	Conductive Adhesive	Yes	300 electrodes/case, 3/bag, 100 bag/case
2284	3-leadwire	.9" x .7"	18" leadwire	Soft Cloth	Conductive Adhesive	Yes	300 electrodes/case, 3/bag, 100 bag/case

To learn more about the 3M™ Red Dot™ Electrodes,  
call the **3M Health Care Helpline** at **1-800-228-3957**  
or visit [www.3M.com/reddotelectrodes](http://www.3M.com/reddotelectrodes).



**3M Health Care**  
**Infection Prevention Division**  
3M Center, Building 275-4E-01  
St. Paul, MN 55144-1000  
USA  
1-800 228 3957  
[www.3M.com/infectionprevention](http://www.3M.com/infectionprevention)

Please recycle. Printed in USA.  
© 3M 2014. All rights reserved.  
70-2010-7538-2 9676HB

3M and Red Dot are trademarks of 3M.  
YokeMate LWS is a registered trademark of  
Advantage Medical Electronics, Inc.