

Hypodermic

Needles & Syringes

Pro Advantage Needles and Syringes are designed to offer these benefits:

Needles

- Sharp, tri-beveled, stainless steel needle with non-coring design
- Unique grinding and polishing process enhances sharpness for ease of penetration
- Smooth surface with light silicone coating to minimize patient discomfort
- Color-coded hubs for easy identification

Syringes

- Latex-free, transparent barrel contains bold graduation markings for accurate filling
- Smooth plunger movement and large, stable finger grips to ensure ease of use
- Double seal stopper to prevent medication leakage
- Backstop to prevent accidental plunger withdrawal



Sterile

Hypodermic Needles


Item #	Description	Hub Color	Box
P929181	18G x 1" Hypodermic Needle	Pink	100
P929185	18G x 1½" Hypodermic Needle	Pink	100
P929201	20G x 1" Hypodermic Needle	Yellow	100
P929211	21G x 1" Hypodermic Needle	Green	100
P929215	21G x 1½" Hypodermic Needle	Green	100
P929221	22G x 1" Hypodermic Needle	Black	100
P929225	22G x 1½" Hypodermic Needle	Black	100
P929231	23G x 1" Hypodermic Needle	Blue	100
P929258	25G x ⅝" Hypodermic Needle	Orange	100
P929251	25G x 1" Hypodermic Needle	Orange	100
P929255	25G x 1½" Hypodermic Needle	Orange	100
P929271	27G x 1" Hypodermic Needle	Grey	100
P929275	27G x 1½" Hypodermic Needle	Grey	100



Luer Lock

Syringes

Item #	Description	Box
P929003	3mL Luer Lock Syringe	100
P929005	5mL Luer Lock Syringe	100
P929010	10mL Luer Lock Syringe	100
P929020	20mL Luer Lock Syringe	100
P929030	30mL Luer Lock Syringe	100
P929050	50mL Luer Lock Syringe	50



Luer Lock

Syringes with Needles

Item #	Description	Box
P923211	3mL Luer Lock Syringe w/ 21G x 1" Needle	100
P923221	3mL Luer Lock Syringe w/ 22G x 1" Needle	100
P923225	3mL Luer Lock Syringe w/ 22G x 1½" Needle	100
P923231	3mL Luer Lock Syringes w/ 23G x 1" Needle	100
P923251	3mL Luer Lock Syringe w/ 25G x 1" Needle	100
P923258	3mL Luer Lock Syringe w/ 25G x ⅝" Needle	100
P921021	10mL Luer Lock Syringe w/ 21G x 1" Needle	100



Tuberculin

Syringes

Item #	Description	Box
P929001	1mL TB Syringe, Luer Slip	100
P921258	1mL Luer Slip TB Syringe w/ 25G x ⅝" Needle	100
P921265	1mL Luer Slip TB Syringe w 26G x ½" Needle	100



For your convenience we also carry other Pro Advantage injection supplies.



Item #	Description	Size	UOM
ALCOHOL PREP PADS			
Sterile			
P902045	Alcohol Prep Pad	1 3/16" x 2 5/8"	200 pk/bx



Item #	Description	Size	UOM
FABRIC ADHESIVE BANDAGES			
Sterile			
P150110	Knuckle Bands	1 1/2" x 3"	100/bx
P150115	Finger Tip	1 3/4" x 2"	100/bx
P150120	Strips	1" x 3"	100/bx
P150125	Strips	2" x 4"	50/bx
P150130	Strips	3/4" x 3"	100/bx



Item #	Description	Size	UOM
SHEER ADHESIVE BANDAGES			
Sterile			
P150135	Strips	1" x 3"	100/bx
P150140	Strips	2" x 4"	50/bx
P150145	Strips	3/4" x 3"	100/bx
P150150	Spots	7/8" Dia.	100/bx



Distributed by: