



NOVUM

Medical Products, Inc.

New Focus. New Commitment. New Innovative Solutions.

Bassinet Model NB-SSxB

Combines beauty with stability. Suited for all-around service in any nursery. Constructed of reinforced, heavy-gauge polished stainless steel for safety and durability.

Features:

- Stainless steel shelf is attached 8" from the ground
- NK Whisper-Glidesm 3" Casters
- Trendelenburg may be created from either end of the bassinet
- Extender shelf option available



Pictured with optional acrylic plastic basket and mattress pad



Dimensions: 31" L x 17.5" W x 37.75" H

Phone | 716.759.7200
Toll Free | 800.274.2742
Fax | 716.759.0700
Email | info@novummed.com
www.novummed.com