



### Training Pads II

Item # 989803139271

Training Pads II can be used with the HeartStart FRx or the HeartStart FRx Trainer on a manikin connected with an external manikin adapter (M5089A) or on a pads placement guide (M5090A for adults; 989803139281 for infants and children). When training on the HeartStart FRx Defibrillator, simply plug in the Training Pads II to activate the FRx's training mode and suspend the defibrillator's ability to deliver a shock.

Users can run eight different training scenarios. Plus, one set of Training Pads II can be used to train for adult and infant/child\* resuscitation. The defibrillator's voice instructions, including CPR Coaching, will be appropriate for treating the simulated victim.

Each training pads case contains a pair of reusable adhesive pads with integrated cable and connector. The case is packaged with an adult pads placement guide, an illustrated guide for using the pads, and Instructions for Use. Replacement Training Pads II (989803139291) are also available, consisting of just a set of pads/wire/plug. An infant/child pads placement guide is also available (989803139281).



### HeartStart FRx Trainer

Item # 861306

This training version of HeartStart FRx Defibrillator includes a Trainer, Training Pads II, a Carry Case, Instructions for Use, Quick Reference, and an external manikin adapter. This product can be used on a manikin equipped with an internal manikin adapter (M5088A) or an external manikin adapter (M5089A).

### Instructor Training Toolkit

Item # 989803139321

DVD only Item # 989803139341

The Training Toolkit includes a DVD and a CD with instructional aids and presentations for teaching groups of people to operate the HeartStart FRx Defibrillator. A DVD training disc can also be purchased separately from the toolkit.



### Infant/Child Key

Item # 989803139311

The Infant/Child Key may be used to simulate infant/child therapy during training with the HeartStart FRx Trainer or with the HeartStart FRx AED with Training Pads II installed.

\* To facilitate training on an infant or child, the Infant/Child Key (989803139311) must be engaged.