



DATA SHEET

HydraFlock® Flocked Swab - Dry Transport System

Mini-tip, polystyrene handle - sterile

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Swab Handle

- 0.098" Dia. (2.489 mm) tapered to 0.049" (1.245 mm)
- 6.0" Length (152.4 mm)

Swab Tip

- 0.135" Dia. (3.429 mm)
- 0.725" Length (18.415 mm)

Qualities

- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Standards

- CE Marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- For nasopharyngeal and urethral cell collection
- Collection of bacterial or micro-organisms
- Forensic evidence collection
- Specimen collection



Ref#	Packaging	
25-3316-H BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube