



**DATA SHEET**

## PurFlock Ultra® Flocked Swab - Dry Transport System

Elongated tip, polystyrene handle with breakpoint at 30mm - sterile

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. PurFlock Ultra® flocked swabs provide excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

### Dimensions

**Overall Length**

- Tube: 6.062" (153.975 mm)
- Swab: 6" (152.4 mm)

**Swab Handle**

- 0.1" Dia. (2.54 mm) tapered to 0.067" (1.702 mm)
- 5.96" Length (151.384 mm)

**Swab Tip**

- 0.219" Dia. (5.563 mm)
- 0.685" Length (17.399 mm)

### Qualities

- Medical grade quality
- Highly absorbent tip
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

### Standards

- CE Marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

### Applications

- Rapid diagnostic testing
- Forensic evidence collection
- Specimen collection
- For collection of bacteria or microorganisms
- DNA testing



Ref#	Packaging
25-3606-U BT	5/100/1 5 boxes per case, 100 tubes per box, 1 swab per tube