

Calgiswab® Calcium Alginate Tipped Applicator

Wood handle



- Calcium alginate tip

- Rigid wood handle

Actual Size

Description

This specialty applicator is tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. It is bonded securely to a rigid wooden handle by an aqueous based adhesive.

Dimensions

6" Overall Length (152.4 mm)

Tip

■ 0.219" Dia. (5.563 mm)

■ 0.641" Length (16.281 mm)

Handle

■ 0.1" Dia. (2.54 mm)

■ 5.812" Length (147.625 mm)

Standards

■ CE marked for medical use

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of sterility and/or compliance available with each shipment upon request

Qualities

■ Medical grade quality

■ Calcium alginate tip contains no fatty acids, which may be inhibitory to many strains of bacteria

Applications

■ Specimen collection for testing and screening

■ Wound care

Ref #	Packaging
806-WA 5M (non-sterile)	5/10/100 (5 boxes per case; 10 bags per box; 100 swabs per bag)
25-806 1WA (sterile)	10/100/1 (10 boxes per case; 100 individually wrapped packages per box; 1 per package)
25-806 2WA (sterile)	10/100/2 (10 boxes per case; 100 individually wrapped packages per box; 2 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com