

DNA-free Cotton Tipped Applicator

Wood shaft



- USP Grade Cotton Tip

- Northern Maine White Birch Shaft

Actual Size



Puritan
Quality since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch shaft by an aqueous based adhesive.

Tested by an independent lab to the most stringent standards (less than 25 pg of human DNA), certified to be DNA-free, DNase-free, and RNase-free.

Dimensions

6" Overall Length (152.4mm)

Tip

- .188" Dia. (4.78mm)

- .625" Length (15.88mm)

Shaft

- .100" Dia. (2.54mm)

- 5.812" Length (147.62mm)

Standards

- CE marked for medical use

- ISO 9001:2008

- ISO 13485:2003

- FDA GMP's

- Certified to be DNA-free & RNase-free

Qualities

- Medical grade quality

- Absorbent tip

- Soft smooth tip

Applications

Forensics : DNA Testing

- Crime scene investigation

- Buccal cell collection

- Specimen collection

Ref #	Packaging
25-806 1WC FDNA <i>Sterile</i>	10/100/1 (10 boxes per case, 100 packages per box, 1 swab per package)

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com