

Polyester Tipped Applicator

Wood handle



- Standard size tip

- Sturdy wood handle

Actual Size

Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Available sterile or non-sterile.

Dimensions

5.937" Overall Length (150.8mm)

Tip

- 0.188" Dia. (4.775mm)
- 0.687" Length (17.45mm)

Shaft

- 0.100" Dia. (2.54mm)
- 5.812" Length (147.625mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of compliance available with each shipment upon request

Qualities

- Medical grade quality
- Contains no fatty acids or latent sugars
- Sturdy wood handle naturally biodegradable

Applications

- Specimen collection
- Patient care
- Surface collection

Ref #	Packaging
806-WD <i>Non-sterile</i>	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-806 1WD <i>Sterile</i>	10/100/1 (10 boxes per case, 100 individually wrapped swabs per box, 1 per package)
25-806 2WD <i>Sterile</i>	10/100/2 (10 boxes per case, 100 individually wrapped swabs per box, 2 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com