

Polyester Tipped Applicator

Wood handle



- Tapered mini-tip

- Sturdy wood handle

Actual Size

Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Available sterile or non-sterile.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.12" Dia. (3.048 mm)
- 0.5" Length (12.7 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.812" Length (147.625 mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of compliance available with each shipment upon request

Qualities

- Medical grade quality
- Contains no fatty acids or latent sugars
- Sturdy wood handle naturally biodegradable
- Tapered mini-tip for hard-to-reach or small areas

Applications

- Specimen collection
- Patient care
- Surface collection

Ref #	Packaging
25-826 2WD Sterile	10/100/2 (10 boxes per case, 100 individually wrapped swabs per box, 2 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com