

# HydraFlock® Flocked Swabs

Polystyrene handle with molded  
break point, standard tip



- Standard flock tip

- Polystyrene handle with molded break point

## Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens. Patents pending: Europe 102513888.4, US 12/849,343.

## Dimensions

6.04" Overall Length (153.416 mm)

### Tip

- 0.200" Dia. (5.08 mm)
- 0.69" Length (17.526 mm)

### Handle

- 0.100 tapered to 0 .066" Dia. (2.54 tapered to 1.676 mm)
- 6.01" Length (152.654 mm)

## Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

## Qualities

- Medical grade quality
- Highly absorbent capacity
- Soft smooth tip

## Applications

- For collection and release of biological specimens
- Buccal cell collection
- Specimen collection
- Forensic evidence collection
- Rapid diagnostic testing

Ref #	Packaging
<b>3306-H</b> <i>Non-sterile</i>	10/100 (10 bags per case, 100 per bag)
<b>25-3306-H</b> <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 swab per package)
<b>25-3306-H TT MC</b> <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped transport devices, 1 swab per device)



**HydraFlock®**

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)  
207-876-3311

Fax 800-323-4153  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
[puritanmedproducts.com](http://puritanmedproducts.com)