


# HydraFlock® Flocked Swabs

Polystyrene handle with molded  
break point, micro ultrafine tip



- Micro ultrafine flock tip

- Polystyrene handle with molded break  
point

Actual Size

## Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens. Patents pending: Europe 102513888.4, US 12/849,343.

## Dimensions

6.03" Overall Length (153.162 mm)

### Tip

■ 0.140" Dia. (3.356 mm)

■ 0.312" Length (7.925 mm)

### Handle

■ 0.100 tapered to 0.036" Dia.  
(2.54 tapered to 0.914 mm)

■ 5.99" Length (152.146 mm)

## Standards

■ CE marked for medical use

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of sterility and/or conformance available  
with each shipment upon request

## Qualities

■ Medical grade quality

■ High absorbent tip

■ Soft smooth tip

## Applications

■ For Nasopharyngeal and Urethral cell collection

■ Rapid diagnostic testing

■ Specimen collection

■ Forensic evidence collection

Ref #	Packaging
<b>3318-H</b> <i>Non-sterile</i>	10/100 (10 bags per case, 100 per bag)
<b>25-3318-H</b> <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 swab per package)
<b>25-3318-H TT MC</b> <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped transport devices, 1 swab per device)



**HydraFlock®**

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)  
207-876-3311

Fax 800-323-4153  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
[puritanmedproducts.com](http://puritanmedproducts.com)