



DATA SHEET

Cotton Tipped Applicator

Doubled-ended cotton tip, one tapered & one regular size, wood handle.

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch shaft by an aqueous based adhesive.

Dimensions

6.0" Overall Length (152.399 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.812" Length (147.625 mm)

Tip

- 0.188" Dia. (4.775 mm) Regular tip
- 0.085" Dia. (2.159 mm) Tapered Tip
- 0.625" Length (15.875 mm)

Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For precision applications: Jewelry Manufacturing and Art Restoration
- General patient care
- Double-ended
- Removing dirt or resins

Ref#	Packaging
821-WC DBL (Non-sterile)	10/5/100 (10 boxes per case, 5 bags per box, 100 per bag)
821-WC (Non-sterile)	10/5/100 (10 boxes per case, 5 bags per box, 100 per bag)

