

# Cotton Tipped Applicator

Wood handle



- USP Grade Cotton, small tip

- Sturdy wood handle

Actual Size

## Description

This cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Small cotton tip ideal for hard-to-reach or small areas.

## Dimensions

■ 5.875" Overall Length (149.225 mm)

### Tip

■ 0.141" Dia. (3.581 mm)

■ 0.437" Length (11.1 mm)

### Handle

■ 0.083" Dia. (2.108 mm)

■ 5.812" Length (147.625 mm)

## Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's

## Qualities

- Medical grade quality
- Small tip for hard-to-reach or small areas
- Sturdy wood handle
- Features minimal cotton for small or hard-to-reach areas

## Applications

**Critical Environment / Industrial / Medical / Veterinary**

- General purpose swab
- Cleaning medical devices
- Removing oil or dirt from machine parts
- Veterinary care
- Cleaning sensors & components

Ref #	Packaging
<b>868-WCS</b> <i>Non-sterile</i>	10/10/100 (10 boxes per case; 10 bags per box; 100 swabs per bag)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)  
207-876-3311

Fax 800-323-4153  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
[puritanmedproducts.com](http://puritanmedproducts.com)