

LifeDop™ 350R (Model # L350R)



Summit Doppler Systems offers a full line of Doppler systems with the sensitivity and performance you depend on. The **LifeDop 350R** large display, table-top Doppler with recharger has excellent sound quality and accepts the full line of non-directional probes for obstetrical or vascular applications. This system is ergonomically designed to comfortably carry and allows easy access to each control feature. Each system is individually tested and inspected to ensure the highest quality standards.

Features:

- Large display for indication of heart-rate and battery status
- Two probe holders for multiple applications
- Superior sound quality from high-fidelity speaker
- Powered by 3 rechargeable batteries
- Extended battery life
- Enhanced probe sensitivity
- Interchangeable probe capability
- User replaceable cable
- Each main unit comes with choice of one probe
- Wall mount available
- Dimensions (h w l): 200 x 150 x 100 mm
- Weight: 1000 grams
- Full 2 year product warranty

Summit Doppler Systems, Inc.
4680 Table Mountain Dr. #150
Golden, Colorado 80403
303-423-7572 or 800-554-5090
Fax # 303-940-7165
www.summitdoppler.com