

SurFlash® and SURFLO® I.V. Catheters



SurFlash® Polyurethane I.V. Catheters

SurFlash offers all the features you value in SURFLO I.V. Catheters with the added benefit of visual confirmation of catheter placement during venipuncture, which occurs in a flash! Other catheter brands only confirm needle placement - *not* catheter placement. SurFlash is made of extremely flexible, kink-resistant and durable polyurethane. With Terumo's unique flashback feature, your staff will attain a higher visual certainty of catheter placement even during difficult venipunctures.

SURFLO® ETFE I.V. Catheters


An ultra-sharp double-bevel introducer needle plus a medical grade lubricant allow for easier penetration and smoother travel through tissue. Preview chamber in clear sure-grip hub gives immediate indication of proper placement. Flexible thinwall catheter with large inner diameter assures good blood flow. All this plus safe, easy-to-handle, tamper-evident packaging make SURFLO® I.V. Catheters a great clinical value.

These products contain no components made of natural rubber latex.

SurFlash: Flashback Feature

With SurFlash, visual confirmation of catheter placement occurs due to a unique flashback feature. Flashback appears between the specially grooved needle and the catheter, providing medical staff with a higher certainty of proper catheter placement. The specially grooved needle also improves the ability to visualize flashback in patients where venipuncture is more difficult. And because SurFlash is made of polyurethane, which is firm at room temperature and softer at body temperatures, patient comfort is enhanced, without compromising proper technique.

NOTE: The needle insertion angle in the illustrations differs from actual insertion angle. When using SurFlash for venipuncture, follow the same procedures and take the same safety precautions as when using Terumo's conventional IV catheter models.



SurFlash Polyurethane I.V. Catheters	
Packed: 50 per box, 4 boxes per case (200).	
Product Code	Description
SR*FF1451	14G x 2" (Orange)
SR*FF1651	16G x 2" (Gray)
SR*FF1664	16G x 2 1/2" (Gray)
SR*FF1832	18G x 1 1/4" (Green)
SR*FF1851	18G x 2" (Green)
SR*FF1864	18G x 2 1/2" (Green)
SR*FF2025	20G x 1" (Pink)
SR*FF2032	20G x 1 1/4" (Pink)
SR*FF2051	20G x 2" (Pink)
SR*FF2225	22G x 1" (Blue)
SR*FF2419	24G x 3/4" (Yellow)

SURFLO ETFE I.V. Catheters	
Packed: 50 per box, 4 boxes per case (200).	
Product Code	Description
SR-OX1451CA	14G x 2" (Orange)
SR-OX1632CA	16G x 1 1/4" (Gray)
SR-OX1651CA	16G x 2" (Gray)
SR-OX1664CA	16G x 2 1/2" (Gray)
SR-OX1832CA	18G x 1 1/4" (Green)
SR-OX1851CA	18G x 2" (Green)
SR-OX1864CA	18G x 2 1/2" (Green)
SR-OX2025CA	20G x 1" (Pink)
SR-OX2032CA	20G x 1 1/4" (Pink)
SR-OX2051CA	20G x 2" (Pink)
SR-OX2225CA	22G x 1" (Blue)
SR-OX2419CA	24G x 3/4" (Yellow)

Call your Terumo representative at 1.800.888.3786 or visit www.terumoTMP.com and click on Find Your Local Terumo Rep.

TerumoTMP.com/IVCatheters

