



*Shown with 044 Accessory

Head Centric Procedure Chair with Programmable Hand Control: Model 4010-650-100

CHAIR SPECIFICATIONS

Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat
Max. 37" at seat (93.9cm)
Seat Width: 28" (71.12cm)
Back Width: 24.5" (62.23cm)
Length: 70"-81" with foot section extended
(177.8cm – 205.74cm)
Leg Rest: 28" W x 13" D (71.12cm W x 33.02cm D)
Dual Electro-Mechanical Actuator System for Efficient Speed
Electrical: 120 VAC, 60HZ
International Model 4010-650-101: 220/240V, 50/60Hz

KEY FEATURES

- Chair Designed for Head and Neck Centric Practices
- 650 lb. (294.84 kg) Patient Weight Capacity
- Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- 4 Power Functions: Hi-Lo, Back, Tilt, Leg Lift & Extension
- Hand Control with OneTouch Patient Positioning®
2 Programmable Positions & Return-to-Home
- Double-articulating Headrest with Height Adjustability and Interchangeable Headrest Options
- Slim Back Cushion for Accessibility
- Removable Seamless Upholstered Cushion
- Trendelenburg Position
- Hospital-Grade Electrical Receptacle with Child Safety Cover (Patient Left)
- Concealed Paper Roll Holder
- 3 Year Limited Warranty

ACCESSORIES

044 - Adjustable/Self Leveling Arms Pair- Field Installable
050 - Removable U-shaped Headrest
051 - Removable Rectangular Headrest
467 - Mounting Plate for Welch Allyn Light Bracket #48955 Patient Right
468 - Mounting Plate for Welch Allyn Light Bracket #48955 Patient Left
830 - Patient Safety Straps
900 - Lifting Bracket Kit
920 - Easy Installation Kit
Fire retardant foam (CAL133) available

COLORS



Hand
Control