

## 6200 Series Mobile Treatment & Supply Cabinets

### Modular Treatment Cabinet with 2 Drawers, 2 Doors and Kydex Work Surface



All-Steel Construction

2 Drawers

2 Doors

Kydex Top

2" Twin Wheel Casters

Available In All 12 UMF Medical Colors

The UMF Medical 6214 Mobile Treatment Cabinet features two drawers, two doors, Kydex work surface, stainless steel push handle and 2" twin wheel carpet casters.

3 year warranty included.

#### Color Chart:



### Standard Features

- 20 gauge steel construction
- Stainless steel push handle
- Kydex top with rail
- 2" (5.1 cm) twin wheel carpet casters
- Heavy duty concealed door hinges
- Full extension roller bearing drawer glides
- Interchangeable drawer fronts
- 12 accent colors available
- 3 year limited warranty

### Specifications

- Width: 25" (63.5cm)
- Height: 34.25" (87.0cm)
- Depth: 18" (45.7cm)

### Options

- 502 – Double door lock
- 505 – Single drawer lock
- 507 – Additional keys (x2)