

Ultra Comfort Stool, Threaded Height Adjustment

Model 6769

All UMF Medical stools are available in 16 best-in-class seamless upholstery colors with Prefixx® Protective Finish to assist with infection prevention and for full exam room cohesion.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Colors



Key features at a glance

- › Ultra Comfort seamless seat
- › 5 leg black plastic base
- › 2" (5.1cm) twin wheel casters
- › Threaded stem for height adjustment
- › 3-year limited warranty

Specifications

Seat diameter: 16" (40.6cm)
 Min. seat height: 20.5" (52.1cm)
 Max. seat height: 28" (71.1cm)
 Net Weight: 14lb (6.4kg)
 Shipping Weight: 19lb (8.6kg)
 Shipping Dimensions: L23xW24xH8
 (58.4cm x 61cm x 20.3cm)

Accessories

150 – Stool glides – set of 5
 111 – Auto-lock occupied casters