

ENT Chair with Foot Operated Pump



Color Chart and Part Number:

8612RI River Rock	8612SN Sand Grey	8612MO Mocha	8612AD Adobe
8612ST Steel Blue	8612MI Midnight Blue	8612SG Sage	8612ON Onyx

Standard Features

- Sturdy steel frame
- Ultra comfort seamless top with antimicrobial guard
- Available in all UMF Medical upholstery options
- 360 degree rotation with swivel lock
- 10 degree elevation from horizontal position
- Manual foot pedal height adjustment
- Manually adjustable backrest
- Manually adjustable footrest
- Hinged arms
- 250 pound (113.4kg) patient weight capacity
- 3 year limited warranty

Specifications

- **Length:** 71" (180.3cm) fully extended
- **Min. Seat Height:** 23" (58.4cm)
- **Max. Seat Height:** 38" (96.5cm)
- **Overall Width:** 31.5" (80.0cm)
- **Seat Width:** 18" (45.72cm)
- **Arm Length:** 23" (58.4cm)

Options

- 031 - Optional backrest with articulating headrest