

### Bariatric and Parent with Child Phlebotomy Chair

#### Standard Features


- Ultra comfort seamless upholstery with antimicrobial guard
- Accommodates bariatric patients and parent with child
- Available in all upholstery color options
- Padded, flip-up, height adjustable arm board
- Heavy duty, all-welded steel frame
- 750 pound (340.2 kg) patient weight capacity
- Adjustable leg levelers
- 3 year limited warranty

#### Options

- 120 - Side tray

#### Specifications

- **Overall dimensions:**  
**Width:** 48-1/2" (123cm)  
**Depth:** 27-3/4" (70cm)  
**Height:** 39-7/8" (101cm)
- **Seat dimensions:**  
**Height:** 21-1/4" (54cm)  
**Width:** 26-5/8" (68cm)
- **Between armboards:** 30" (76cm)



Shown with optional -120 side tray

#### Color Chart and Part Number:

<b>8690RI</b> River Rock	<b>8690SN</b> Sand Grey	<b>8690MO</b> Mocha	<b>8690AD</b> Adobe
<b>8690ST</b> Steel Blue	<b>8690MI</b> Midnight Blue	<b>8690SG</b> Sage	<b>8690ON</b> Onyx