

### Spacious Storage, Angled Top, Full View Doors



All-Welded Construction  
Hinged, Glass Panel Doors with Lock  
5 Adjustable Glass Shelves  
Adjustable Leveling Glides  
3 Year Limited Warranty

The UMF Medical SS7834 Large stainless steel storage and supply cabinet features an angled top and full height stainless steel hinged doors with glass panes. Locking T-handle on right side door locks both doors in the closed position. Stationary T-handle on left side door. Five removable glass shelves, adjustable in ½" (1.3cm) increments.

The UMF Medical SS7834 Large Stainless Steel Storage and Supply Cabinet comes with a 3 year limited warranty.

#### Standard Features

- Heavy duty, all-welded stainless steel construction
- Five removable, adjustable stainless steel shelves
- Full Height, stainless steel sliding doors with glass panes
- Convenient T-handle door latches with lock
- Adjustable leveling glides
- 3 year limited warranty

#### Specifications

- Width: 35.875" (91.1cm)
- Height: 77.25" (196.2cm)
- Depth: 16" (40.6cm)

#### Options

- 507 – Additional keys (x2)
- NPN – solid panel doors instead of glass panes