

Blood-Drawing Chair 2570/2571



ISO 9001-2008 Certified

Value and quality at an affordable price

Economical, yet designed with durability and longevity in mind. The high quality features of our Blood Drawing Chairs will provide you with years of service.



KEY FEATURES

- Durable plastic seat/back
- Laminated cabinet (2570 only)
- Drawer roller glides (2570 only)
- Sturdy leveling feet
- Adjustable height armrests

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

BLOOD-DRAWING CHAIR 2570/2571

- 1 Durable plastic resin seat endures years of use.
- 2 Flip-up, contoured urethane armrest provides easy access and added safety for the user.
- 3 Adjustable height armrest.
- 4 High impact leveling feet ensure a sturdy chair on almost any surface.
- 5 Fully laminated cabinet and storage drawer.
- 6 Powder-coated steel frame.



SPECIFICATIONS

2570/2571

Seat Height	20.5" (52cm)
Seat Width	18" (45.7cm)
Seat Depth	15" (38cm)
Overall Height	36" (91.4cm)
Overall Width	2570: 35" (88.9cm)
	2571: 29" (73.7cm)
Back Height Above Seat	16" (40.6cm)
Overall Length	20.5" (52cm)
Adjustable Arm Height Range	27"-33" (68.6 - 83.8cm)
Weight	2570: 46 lbs. (20.9kg)
	2571: 32 lbs. (14.5kg)
Weight Capacity	300 lbs. (136kg)

*Weight shown as a minimum-maximum range & varies depending on chair options ordered.

OPTIONS



Pivot Arm **PR**
(both models)
Dual Pivot Arm **P2**
(2571 only)



IV Pole **IV**
(right-rear only)