



*The Winco Phlebotomy Collection*

# Blood-Drawing Chair



Winco's Blood-Drawing Chairs combine functional design with quality and durability. The collection offers functionality, comfort and convenience for both your patients and your staff.

- All Winco blood-drawing chairs feature an adjustable-height, flip-up, steel reinforced, contoured urethane armrest for optimal comfort, easy access, and added safety.
- High-impact, leveling feet provide stable seating on almost any surface.
- Manufactured with sturdy, powder-coated steel frames for maximum strength and years of trouble-free service.

This collection offers an array of features and options to meet your application and patient treatment needs.

### **BASIC BLOOD-DRAWING CHAIR (2570/2571)**

Economical and extremely sturdy; Winco's Basic Blood-Drawing Chairs feature a durable plastic resin seat, adjustable armrests, and non-marring leveling feet. A clean and simple design that stands up to years of use.



### **XL BLOOD-DRAWING CHAIR (2574/2575)**

In addition to adjustable armrests and high impact leveling feet our XL Blood-Drawing Chairs feature a 24-inch wide padded seat, adjustable armrests and are designed with added weight capacity to safely accommodate patients up to 450 lbs.

### **VERSATILE STORAGE CABINET**

High impact aluminum & plastic storage cabinet with a slotted drawer and a contoured urethane armrest are included on blood-drawing chair models 2570, 2572, and 2574 in our Phlebotomy Collection.



### **DESIGNER BLOOD-DRAWING CHAIR (2578)**

Winco's Designer Blood-Drawing Chair can easily be adjusted up or down based on the height requirements of the phlebotomist, doctor, nurse or clinician. Additionally, our Designer Blood-Drawing Chair also features a footrest with a no-slip grip and is available with an optional side tray.

*\* Shown with optional side table*



The Winco Phlebotomy Collection

# Blood-Drawing Chair

## ANATOMY OF THE BLOOD-DRAWING CHAIR\*

Steel reinforced, flip-up, contoured urethane armrest for easy access and added safety

Comfortable seat & back rest

High impact plastic storage drawer features a spill-containing top and a rugged aluminum box. (2570, 2572, & 2574 only)\*\*



Adjustable height armrests

High-impact leveling feet = a sturdy chair on almost any surface

Durable, powder-coated steel frame

Optional pivot arm available on all models\*\*\*

Padded models available in over 40 vinyl colors, visit our website or contact us for representative examples.

\* Shown: Model 2575

\*\* Cabinet and storage drawer available on models 2570, 2572 and 2574.

\*\*\* Optional side table shown

## MODELS & SPECIFICATIONS



Model 2570 shown



Model 2572 shown



Model 2574 shown



Model 2578 shown with optional side table

MODEL	2570 / 2571	2572 / 2573	2574 / 2575	2578
Seat Height	20.5" (52.07 cm.)	22.5" (57.15 cm.)	22.5" (57.15 cm.)	25.5-31" (64.8 - 78.7 cm.)
Seat Width	18" (45.72 cm.)	18.75" (47.63 cm.)	24.25" (61.60 cm.)	19" (48.26 cm.)
Seat Depth	15" (38.10 cm.)	18.5" (46.99 cm.)	18.5" (46.99 cm.)	19" (48.26 cm.)
Overall Height	36" (91.44 cm.)	44.5" (113.03 cm.)	36.5" (92.71 cm.)	47-52.5" (119.38 - 133.4 cm.)
Overall Width*	35.25" (89.5 cm.) / 29" (73.66 cm.)	35.25" (89.5 cm.) / 29" (73.66 cm.)	40.25" (102.2 cm.) / 34" (86.36 cm.)	27.5" (69.85 cm.)
Back Height Above Seat	16" (40.6 cm.)	24.5" (62.2 cm.)	16" (40.6 cm.)	24" (61 cm.)
Overall Length	20.5" (52.07 cm.)	27" (68.58 cm.)	25" (63.50 cm.) / 25" (63.50 cm.)	27" (68.58 cm.)
Adjustable Arm Height Range	27"-33" (68.6 - 83.8 cm.)	28.25"-34.25" (71.76 - 87 cm.)	28.25"-34.25" (71.76 - 87 cm.)	33-44.5" (83.8 - 113 cm.)
Weight	35 lbs. (15.9kg.) / 32 lbs. (14.51 kg.)	56 lbs. (25.4 kg.) / 50 lbs. (22.68 kg.)	61 lbs. (27.7 kg.) / 55 lbs. (24.95 kg.)	**65-73 lbs. (29.5 - 33.11 kg.)
Weight Capacity	300 lbs. (136.08 kg.)	300 lbs. (136.08 kg.)	450 lbs. (204.12 kg.)	300 lbs. (136.08 kg.)

\* Add 10.5" to chair width for each raised side table.

\*\* Weight shown as a minimum-maximum range & varies depending on chair options ordered.

## FEATURES

MODEL	TB-133	SEAT & BACK			CABINET & DRAWER	SIDE TABLE (R)	FOOTREST BAR	PIVOT ARM		
		PLASTIC RESIN	PADDED	HEIGHT ADJ SEAT				SINGLE (PR)	DUAL (P2)	IV POLE
2570	N/A	S	N/A	N/A	S	N/A	N/A	0	N/A	0 (right rear)
2571	N/A	S	N/A	N/A	N/A	N/A	N/A	0	0	0 (right rear)
2572	0	N/A	S	N/A	S	N/A	N/A	0	N/A	0
2573	0	N/A	S	N/A	N/A	N/A	N/A	0	0	0
2574	0	N/A	S	N/A	S	N/A	N/A	0	N/A	0 (right rear)
2575	0	N/A	S	N/A	N/A	N/A	N/A	0	0	0 (right rear)
2578	0	N/A	S	S	N/A	0	S	0	0	N/A

S-Standard 0-Optimal N/A-Not Available

**WARRANTY** Contact Winco or visit our website for warranty details.

**VINYL** Padded models available in over 40 vinyl colors, visit our website or contact us for swatches.

**NOTE** Winco reserves the right to change specifications without notice.



## Contact Us

To speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com



Made in USA