

Privacy Screen

3630/3640/3650/3400 Standard White
3639/3649/3659/3409 Sure-Chek®



ISO 9001-2008 Certified

A privacy screen for every purpose.

Beautiful style, great colors! Winco has a screen to meet your every need; whether it is a light-weight aluminum frame or a space saving wall mounted screen.



3400

KEY FEATURES

- Replaceable screen panels
- Bacteriostatic vinyl

3400 Style

- Pivots 180 degrees
- Extends from 40" - 92.5"

3630 Style

- Hinged screen panels
- Wide base for stability



3639

4 panel (3640/3649) & 5 panel (3650/3659) screens are also available.

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

STEEL FRAME PRIVACY SCREEN

3630/3640/3650/3400 STANDARD WHITE

3639/3649/3659/3409 SURE-CHEK®

- 1 Choose from our bacteriostatic vinyl panels or our OPTIONAL Sure-Chek® antimicrobial material.
- 2 Hinged panel ensures years of problem-free use.
- 3 Wide base for added stability



SPECIFICATIONS

	3630/3639/3640/3649/3650/3659	3400/3409
Overall Height.....	69" (175.3cm)	72" (182.9cm)
Overall Width.....	3630/3639: 60" (152.4cm) 3640/3649: 80" (203.2cm) 3650/3659: 100" (254cm)	40" - 92.5" (101.6 - 235cm)
Single Panel Length.....	46" (116.84cm)	69" (175.26cm)
Single Panel Width.....	20" (50.8cm)	N/A
Weight.....	3630/3639: 9 lbs. (4.1kg) 3640/3649: 12 lbs. (5.4kg) 3650/3659: 15 lbs. (6.8kg)	.8 lbs. (3.6kg)

COLORS



Optional Optional Optional Optional Optional

Printed colors shown may vary slightly from actual color of product.
Sure-Chek® vinyls are antimicrobial, antistatic, odor resistant, stain resistant, & fire retardent.