

Treatment Table w/ face cutout 8510



ISO 9001-2008 Certified

Customize your treatment table!

Our tables offer versatility, comfort and strength. Build onto our basic models using shelves, cabinets, and drawers to meet your specific needs.



8510



8510-SH



8510-C1



8510-D1

KEY FEATURES

- Face cutout w/ paper holder
- Durable foam construction
- Sturdy steel leg frames
- Paper holder
- Adjustable leveling feet

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

TREATMENT TABLE w/ FACE CUTOUT 8510

- 1 Heavy gauge expanded vinyl & medium density foam for superior comfort.
- 2 Standard paper holder.
- 3 High strength powder-coated steel frame.
- 4 High impact leveling feet ensure a sturdy table on almost any hard surface.
- 5 Face cutout with face paper holder and cutter.



SPECIFICATIONS

8510

Overall Height* 30" (76cm) STANDARD*
 Overall Width 28" (71cm)
 Overall Length 72" (183cm)
 Weight** 74 lbs. (33.6kg)**
 Weight Capacity 400 lbs. (181kg)

*Special height legs available in 19", 25", 36"
 **Weight varies depending on table options ordered.
 Add 25 lbs. for 8510-SH

OPTIONS



TB-133 **TB**



*Special Height Legs



Shelf **SH**



Cabinet **C1**



Drawer **D1**



Matching Pillow
MPOO

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS