

UpScale xM450 SPECIFICATIONS

Weight of table	250lb/113kg
Dimensions of shipped table (boxed, flat)	60" x 36" x 36"
Weight of shipped table	320 lb/145kg

Dimensions (as installed)

Upholstery pad dimensions	28" x 55"
Upholstery pad thickness	2.75" (more at head)
Leg support dimensions	18" x 18"
Length of table with scale display arm installed	63" long
Length of table with leg support and scale display arm	80.5" long
Width of table with grab bars	33.5" wide
Height of table – floor to seat	17" minimum 33.5" maximum
Seat back range	fully powered, 0° to 61°
Pelvic tilt range	manual, 0° to 9°
Foot stirrups	manually extended 23" through 3 lateral positions
Knee crutches	manually extended 12.5" high from upholstered seat through infinite positions
Tape measure	120"/305cm. long retractable soft tape
Power cord	96" long

Safe Working Loads

Exam deck (including head, seat and legs)	450 lb/204 kg
Leg support	125 lb/57 kg

Hand Pendant

Fully powered, controls table height and seat back	
Weight	.7 lb/.32 kg

Scale Display	liquid display
Dimensions of unit	6.5" x 3" x 3.5"
Power	4 AA batteries
Accuracy	+/- 0.5%
Paper Roll	21" roll or smaller

Finishes

Frame Body	powder-coated paint over steel
Grab Bars	chrome-polished steel
Knee Crutch	polyurethane foam padding
Upholstery*	vinyl with PreFixx, meets CAL 117
<i>Optional upholstery that meets CAL 133 available as an upgrade</i>	

Electrical

Power Supply	120 VAC, 60 hz
Maximum Draw	2.7 A
Protection against electric shock	Class 1
Classification of patient-applied parts	Type B
Protection against water ingress	
For electrical enclosures	IPX0
Fuse rating	6 A 250 V
Duty cycle (motor run time)	10% Max 2 min on/ /18 min off
Certifications	UL 60601-1; IEC 60601-1; CAN-CSA C22.2 No. 60601-1;
Operating conditions	
Temperature range	41°F to 104°F (5°C to 40°C)
Humidity	20% to 90% @ 86°F/30°C
Transportation/Storage Conditions	
Temperature range	41°F to 104°F (5°C to 40°C)
Relative humidity	10% to 90% (non-condensing)
Atmospheric Pressure	80-106 kPa

*Fire Code Ratings: All upholstery complies with California Bureau of Home Furnishing Technical Bulletin 117 and California Code of Regulations, Sec. 93120-93120.12, Title 17. Optional upholstery is available that complies with California Bureau of Home Furnishing Technical Bulletin 133.