



CircuFlow™ 4-Chamber Garment Measurement Chart

Lower Extremity		Model	Length (in.)		Circumference (in.)		
Size	Type	4 chamber	Inseam	Ankle	Knee	Thigh/Groin**	
Small	Full Leg	D-300S	23	21	25	30	
Small Wide	Full Leg	D-300SW	26	24	29	35	
Small Extra Wide	Full Leg	D-300SXW	24	28	35	40	
Medium	Full Leg	D-300M	29	22	27	31	
Medium Wide	Full Leg	D-300MW	29	25	31	35	
Medium Extra Wide	Full Leg	D-300MXW	28	29	36	40	
Large	Full Leg	D-300L	34	21	27	31	
Large Wide	Full Leg	D-300LW	36	22	31	40	
Half	Half Leg	D-301H	16.5	22	28	-	
Half Wide	Half Leg	D-301HW	19	24	30	-	

Upper Extremity		Model	Length* (in.)		Circumference (in.)		
Size	Type	4 chamber	Axilla to Fingertips	Axilla	Elbow	Wrist	
Medium	Arm	D-302M	26.5	22	19	13	
Large	Arm	D-302L	32	23	20	14	

Lower Extremity		Model	Length (cm.)		Circumference (cm.)		
Size	Type	4 chamber	Inseam	Ankle	Knee	Thigh/Groin**	
Small	Full Leg	D-300S	58.4	53.3	63.5	76.2	
Small Wide	Full Leg	D-300SW	66	61	73.7	88.9	
Small Extra Wide	Full Leg	D-300SXW	61	71.1	88.9	101.6	
Medium	Full Leg	D-300M	73.7	55.9	68.6	78.8	
Medium Wide	Full Leg	D-300MW	73.7	63.5	78.7	88.9	
Medium Extra Wide	Full Leg	D-300MXW	71.1	73.7	91.4	101.6	
Large	Full Leg	D-300L	86.4	53.3	68.6	78.7	
Large Wide	Full Leg	D-300LW	91.4	55.9	78.7	101.6	
Half	Half Leg	D-301H	41.9	55.9	71.1	-	
Half Wide	Half Leg	D-301HW	48.3	61	76.2	-	

Upper Extremity		Model	Length* (cm.)		Circumference (cm.)		
Size	Type	4 chamber	Axilla to Fingertips	Axilla	Elbow	Wrist	
Medium	Arm	D-302M	67.3	55.9	48.3	33	
Large	Arm	D-302L	81.3	58.4	50.8	35.6	

* Length corresponds to the entire sleeve length from axilla to fingertips ** The circumference of the upper rim of the sleeve



CircuFlow™ 8-Chamber Garment Measurement Chart

Lower Extremity		Model	Length (in.)		Circumference (in.)		
Size	Type	8 chamber	Inseam	Ankle	Knee	Thigh/Groin**	
Small	Full Leg	D-380S	23	21	25	30	
Small Wide	Full Leg	D-380SW	26	24	29	35	
Small Extra Wide	Full Leg	D-380SXW	24	28	35	40	
Medium	Full Leg	D-380M	29	22	27	31	
Medium Wide	Full Leg	D-380MW	29	25	31	35	
Medium Extra Wide	Full Leg	D-380MXW	28	29	36	40	
Large	Full Leg	D-380L	34	21	27	31	
Large Wide	Full Leg	D-380LW	36	22	31	40	

Upper Extremity		Model	Length* (in.)		Circumference (in.)		
Size	Type	8 chamber	Axilla to Fingertips	Axilla	Elbow	Wrist	
Medium	Arm	D-382M	26.5	22	19	13	
Large	Arm	D-382L	32	23	20	14	
Medium	Arm & Shoulder	D-383M	26.5	22	19	13	
Large	Arm & Shoulder	D-383L	32	23	20	14	

Lower Extremity		Model	Length (cm.)		Circumference (cm.)		
Size	Type	8 chamber	Inseam	Ankle	Knee	Thigh/Groin**	
Small	Full Leg	D-380S	58.4	53.3	63.5	76.2	
Small Wide	Full Leg	D-380SW	66	61	73.7	88.9	
Small Extra Wide	Full Leg	D-380SXW	61	71.1	88.9	101.6	
Medium	Full Leg	D-380M	73.7	55.9	68.6	78.8	
Medium Wide	Full Leg	D-380MW	73.7	63.5	78.7	88.9	
Medium Extra Wide	Full Leg	D-380MXW	71.1	73.7	91.4	101.6	
Large	Full Leg	D-380L	86.4	53.3	68.6	78.7	
Large Wide	Full Leg	D-380LW	91.4	55.9	78.7	101.6	

Upper Extremity		Model	Length* (cm.)		Circumference (cm.)		
Size	Type	8 chamber	Axilla to Fingertips	Axilla	Elbow	Wrist	
Medium	Arm	D-382M	67.3	55.9	48.3	33	
Large	Arm	D-382L	81.3	58.4	50.8	35.6	
Medium	Arm & Shoulder	D-383M	67.3	55.9	48.3	33	
Large	Arm & Shoulder	D-383L	81.3	58.4	50.8	35.6	

* Length corresponds to the entire sleeve length from axilla to fingertips ** The circumference of the upper rim of the sleeve